

# EndoTNeedle™

Endoscopic Injection Needle

 **GI Supply®**  
Specialty Endoscopic Products



Click It & Stick It



# A tough needle, at a fair price

## Your new everyday needle

EndoTNeedle was engineered to balance performance and value while delivering a highly reliable needle for your most tortuous cases. Designed with quality materials, features from premium needles, and a handle your staff will love. The EndoTNeedle is your new everyday needle.



**Single hand  
extension**




**Click to  
retract**

## Click It

- Single-handed use in an ergonomic handle
- Thumb activated needle extension and spring loaded retraction



**Reliably works in  
retroflex**



**A sharp needle  
and soft hub**

## **Stick It**

- Reliable performance in retroflex and tortuous paths
- Stainless steel sheath prevents kinking
- Smooth round hub which aids in the closure and blockage needed to stop GI bleed cases
- Ultra sharp needle





## The perfect fit, in a kit!

Spot® Ex and EndoTNeedle go together perfectly. Now you can order both with one item number. The kit includes a box of Spot Ex and a box of the EndoTNeedle.

Item number GIS-4501.

Add the EndoTNeedle on your next GI Supply order!

Order Number	Minimum Working Channel	Working Length	Needle Extension	Needle Gauge	Units/Box
ET1300	2.8mm	160cm	4mm	25	10
ET1301	2.8mm	230cm	4mm	25	10
ET1302	2.8mm	160cm	5mm	25	10
ET1303	2.8mm	230cm	5mm	25	10
ET1320	2.8mm	160cm	4mm	23	10
ET1321	2.8mm	230cm	4mm	23	10
ET1322	2.8mm	160cm	5mm	23	10
ET1323	2.8mm	230cm	5mm	23	10